

LASER SPECTRUM ANALYZER

772 Series




Spectral analysis of pulsed IR and Mid-IR lasers.

The 772B-MIR Laser Spectrum Analyzer from Bristol Instruments combines proven Michelson interferometer technology with fast Fourier transform analysis to characterize the spectral properties of lasers that operate from 1 to 12 µm. What makes this system so unique is that it employs a sophisticated algorithm to enable the measurement of pulsed lasers that have a repetition rate as low as 100 Hz.

The model 772B-MIR collects laser pulses over time to "build" an interferogram that is sufficient to convert to a spectrum. This spectrum has a resolution of 4 GHz and wavelength accuracy of \pm 10 parts per million (\pm 0.08 nm at 8 μ m). An optical rejection ratio of greater than 20 dB is achieved assuming a sufficient number of pulses (~30,000) are used to generate the interferogram.

Pulse Rate: 10 kHz Pulse Length: 50 ns g 10




Spectra of QCL operating with a pulse length of 50 ns and at various repetition rates.

Key Features:

- Spectral analysis of pulsed and CW lasers.
- Operation from 1 to 12 µm.
- Spectral resolution of 4 GHz.
- Wavelength accuracy as high as ± 0.01 nm.
- Optical rejection ratio as high as 20 dB.

PRELIMINARY SPECIFICATIONS



SPECIFICATIONS	772 Series
MODEL	772B-MIR
LASER TYPE ¹	CW and pulsed (repetition rate > 100 Hz)
WAVELENGTH	
Range	1 - 12 µm
Accuracy ^{2, 3.}	± 10 ppm
	\pm 0.08 nm @ 8 μm \pm 0.0125 cm $^{-1}$ @ 1250 cm $^{-1}$ \pm 875 MHz @ 37,500 GHz
Spectral Resolution ⁴	4 GHz
Calibration	Continuous - built-in standard HeNe laser
Display Resolution	8 digits
Units ⁵	nm, μm, cm·¹, GHz, THz
OPTICAL REJECTION RATIO ⁶	> 20 dB (> 30,000 pulses)
MINIMUM INPUT POWER 6	0.005 - 2.5 μW
MEASUREMENT TIME 6	Approximately 2x time required to collect chosen number of pulses (i.e., 6 s for 30,000 pulses from 10 kHz repetition rate laser)
INPUTS/OUTPUTS	
Optical Input ⁷	Collimated beam 3 mm diameter aperture, visible tracer beam to facilitate alignment
Instrument Interface	USB and Ethernet with Windows-based display program Library of commands (SCPI) for custom and LabVIEW programming using any PC operating system
COMPUTER REQUIREMENTS 8	PC running Windows 10, 1 GB available RAM, USB 2.0 (or later) port, monitor, pointing device
ENVIRONMENTAL ⁶	
Warm-Up Time	None
Temperature Pressure Humidity	$+15^{\circ}$ C to $+30^{\circ}$ C (-10° C to $+70^{\circ}$ C storage) 500 − 900 mm Hg ≤ 90% R.H. at $+40^{\circ}$ C (no condensation)
DIMENSIONS AND WEIGHT	
Dimensions (H x W x D)	7.5" x 6.5" x 15.0" (191 mm x 165 mm x 381 mm)
Weight	14 lbs (6.3 kg)
POWER REQUIREMENTS	90 - 264 VAC, 47 - 63 Hz, 50 VA max
WARRANTY	5 Years (parts and labor)

- (1) Specifications are for pulsed laser mode of operation. For CW lasers, operation and specifications are identical to model 771B-MIR.

- 2) Defined as measurement uncertainty, or maximum wavelength error, with a confidence level of ≥ 99.7%.
 3) Wavelength axis is calibrated to system's accuracy specification.
 4) Defined as the measured full width at half maximum intensity (FWHM) of an infinitely narrow optical signal.
- (5) Data in units of nm, µm, and cm⁻¹ are given as vacuum values.
 (6) Characteristic performance, but non-warranted.
- (7) MIR required beam height is 5.4 ± 0.25 ".
- (9) For use with Windows-based display program. Interface with SCPI can be done using any PC operating system.
 (9) MIR instrument height is adjustable (7.25 ± 0.25") for alignment purposes.

 $Bristol\ Instruments\ reserves\ the\ right\ to\ change\ the\ specifications\ as\ may\ be\ required\ to\ permit\ improvements\ in\ the\ design\ of\ its\ products.$ Specifications are subject to change without notice.



