


Beamsplitters



Beamsplitters Selection Guide


B045

Half Mirror Cube



Non-polarizing Cube Half Mirrors
NPCH

B046



Hybrid Cube Half Mirrors
HBCH

B048



Chromium Cube Half Mirrors
CSCH

B049



Dielectric Cube Half Mirrors
CSMH

B050

Half Mirror Plate



Ultra Broadband Dielectric Half Mirrors
PMH/PSMH

B052



Thin Plate Beamsplitter
MPSMH

B054



Laser Line Plate Half Mirrors
PSMH

B055



Chromium Plate Half Mirrors
PSCH

B058

Application Note

B059

Beamsplitters




Dielectric Cube Beamsplitters
CSM

B060



Dielectric Plate Beamsplitters
PSM


B061



Variable Beamsplitter
Light path corrector
VBS/WSQNA/WBNA

B062

Harmonic Separators



Harmonic Separators
YHS

B064

Beam Samplers



Beam Samplers
BS4

B066

Others



Polka Dot Beamsplitter
PDBS

B067



Contact sheet for Special Half Mirror
and Beamsplitter and Beam Sampler


B068

Beamsplitters selection Guide

A beamsplitter is an optic that splits light into 2 directions. The split ratio of light transmittance and reflectance is 1:1 and is called a half mirror.








The 2 forms of beamsplitters are cube and plate type.

Type	Overview drawing	Features	Application
Plate		Good fit for large beam size applications at a reasonable price. Advantages are: minimal back reflection, compact light-path as compared to cube type beamsplitters and low chromatic dispersion. There may be a slight offset of the transmitted beam due to refraction. For 45 degrees incident application, the clear aperture would be elliptical. There may be some vignetting on angle of incidence.	Large beam size optical set up. Used in large beam size optical layouts. Used for monitoring optical systems, split beams into different wavelengths, polarizations or intensities.
Cube		Can be applied at its maximum effective area from any incident direction, easy to be applied in optical design and simple for optical set up adjustment. High cost and high weight for large beam size application. Feedback light at less than 1% may happen. The transmittance light through the cube is longer than a plate type, the chromatic dispersion is higher. Eliminates the problem of beam deviation.	For a compact size optical set up. For high accuracy experiment and optical set up usage.



Experimentation with laser (Linear polarized light)

Lasers are used to evaluate our half mirrors and with the polarization properties of the laser, we are able to check the change of light splitting ratios.


Type	Image	Application	Polarization constraints received when laser light is used	Polarization dependency
Non-polarizing (NPCH) Reference B046		For high accuracy laser experiment with accurate light ratios at any polarization levels.	The light ratios at 1:1 stay stable even when the polarization situation changes. No power loss.	Small
Hybrid (HBCH) Reference B048		For multi-wavelength light splitting solutions.	Light ratio at 1:1 from any specified light incident direction will remain similar.	
Laser Line Plate (PSMH) Reference B055		Large beam size, multi mirror optical set up with small power light source and supports high power laser light splitting.	Polarization at 45 degree (AOI) or circle polarization light with no power loss detected.	
Chromium Plate (PSCH) Reference B058		Large beam size and observation optical system.	Polarization at 45 degree (AOI) or circle polarization light with 36% absorption of light power.	
Chromium Cube (CSCCH) Reference B049		For basic laser experiments and compact optical solutions. Great entry level price.	Polarization at 45 degree (AOI) or circle polarization light with 40% absorption of light power.	
Dielectric Cube (CSMH) Reference B050		For general white light and non-polarizing light i.e. LED light splitting solutions.	Polarization at 45 degree (AOI) or circle polarization light with no power loss detected.	

Non-polarizing Cube Half Mirrors | NPCH


RoHS

Half mirrors have an even 1:1 ratio of reflection and transmission in both linear polarized light and normal light source.

- The reflection to transmission ratio is 1:1 regardless of the polarization condition from the input beam.
- Depending on polarization, the reflection to transmission ratio of these products does not vary.
- The laser lines listed correspond to various laser wavelengths.
- Narrowband multi-layer AR coatings are applied to the four surfaces of the cube.
- These cubes are designed for a single wavelength because the effective bandwidth of a non-polarizing coating is narrow.




Schematic



Outline Drawing

(in mm)



- Tolerance
- A ± 0.2
- B ± 0.2
- C ± 0.1

Specifications

Material	BK7, Synthetic fused silica
Surface flatness of substrate	$\lambda/4$
Beam Deviation	<5'
Coating	Hypotenuse Surface: Dielectric multi-layer nonpolarizing coating Four Surfaces: Multi-layer anti-reflection coating
Incident angle	0°
Divergence ratio (reflectance : transmittance)	1 : 1
Laser Damage Threshold	0.3J/cm ² (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	20-10
Clear aperture	85% of Circle to actual dimension (80% of actual aperture for 5 and 7mm dimension (A=B=C) products.)

Guide

- ▶ Please contact our Sales Team for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ Non-polarizing beam splitter (plate type) is also possible.
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division.

Attention

- ▶ Input beam from the prism side is indicated by a "O".
- ▶ Phase retardation of light input will not be preserved. Use a wave-plate for phase compensation.
- ▶ Wavelength dispersion of transmitted and reflected light is derived from refractive index and glass thickness and when diverging or introducing a focusing beam, chromatic aberration or spherical aberration may occur.

266 – 532nm

Part Number	Wavelength Range [nm]	A=B=C [mm]	Material	Transmittance [%]
NPCH-10-2660	266	10	Synthetic fused silica	50±10
NPCH-15-2660	266	15	Synthetic fused silica	50±10
NPCH-20-2660	266	20	Synthetic fused silica	50±10
NPCH-10-3550	355	10	Synthetic fused silica	50±7
NPCH-15-3550	355	15	Synthetic fused silica	50±7
NPCH-20-3550	355	20	Synthetic fused silica	50±7
NPCH-10-4050	405	10	BK7	50±7
NPCH-15-4050	405	15	BK7	50±7
NPCH-20-4050	405	20	BK7	50±7
NPCH-10-4880	488	10	BK7	50±5
NPCH-15-4880	488	15	BK7	50±5
NPCH-20-4880	488	20	BK7	50±5
NPCH-10-5145	514.5	10	BK7	50±5
NPCH-15-5145	514.5	15	BK7	50±5
NPCH-20-5145	514.5	20	BK7	50±5
NPCH-10-5320	532	10	BK7	50±5
NPCH-15-5320	532	15	BK7	50±5
NPCH-20-5320	532	20	BK7	50±5

Compatible Optic Mounts

PLH-25, -40 / KKD-25PHRO, -40PHRO

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others




Cube Type: Nonpolarizing 632.8 – 1550nm

Part Number	Wavelength Range [nm]	A=B=C [mm]	Material	Transmittance [%]
NPCH-05-6328	632.8	5	BK7	50±5
NPCH-10-6328	632.8	10	BK7	50±5
NPCH-15-6328	632.8	15	BK7	50±5
NPCH-20-6328	632.8	20	BK7	50±5
NPCH-10-6700	670	10	BK7	50±5
NPCH-15-6700	670	15	BK7	50±5
NPCH-20-6700	670	20	BK7	50±5
NPCH-10-7800	780	10	BK7	50±5
NPCH-15-7800	780	15	BK7	50±5
NPCH-20-7800	780	20	BK7	50±5
NPCH-10-8300	830	10	BK7	50±5
NPCH-15-8300	830	15	BK7	50±5
NPCH-20-8300	830	20	BK7	50±5
NPCH-10-10640	1064	10	BK7	50±5
NPCH-15-10640	1064	15	BK7	50±5
NPCH-20-10640	1064	20	BK7	50±5
NPCH-10-13000	1300	10	BK7	50±5
NPCH-15-13000	1300	15	BK7	50±5
NPCH-20-13000	1300	20	BK7	50±5
NPCH-10-15500	1550	10	BK7	50±5
NPCH-15-15500	1550	15	BK7	50±5
NPCH-20-15500	1550	20	BK7	50±5


Typical Transmittance Data

T: Transmission


NPCH-2660




NPCH-3550




NPCH-4050




NPCH-4880




NPCH-5145




NPCH-5320




NPCH-6328




NPCH-6700




NPCH-7800




NPCH-8300




NPCH-10640



NPCH-13000



NPCH-15500



Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers


Others

Hybrid Cube Half Mirrors | HBCH


RoHS Catalog Code W3013

Low polarizing cube half mirrors can be used for broadband visible and infrared light. Cubes are applicable for use in polarizing systems and lasers with multiple wavelength or visible light.


- This hybrid coating is a special combination of metallic and dielectric multi-layers that result in low polarizing and broadband coating.
- As it is cube shaped, there will not be any lateral shift of the optical axis when a normal incident beam is applied. During transmission and reflection of lights, the aperture remains unchanged.
- Even when the orientation of linear polarization has been changed, beams are equally divided as reflected (R) : transmitted (T) (ratio is 1:1)



Schematic



Outline Drawing




Specifications	
Material	BK7
Surface flatness of substrate	$\lambda/4$
Beam Deviation	<5'
Coating	Hypotenuse surface: Hybrid coating (dielectric multi-layer coating and metallic coating) Four surfaces: Multi-layer anti-reflection coating
Incident angle	0°
Divergence ratio (reflectance : transmittance)	1 : 1
Laser Damage Threshold	0.3J/cm ² (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	40-20
Clear aperture	85% of actual dimension

- ### Guide
- ▶ Please contact our Sales Team for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
 - ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division.

- ### Attention
- ▶ Input beam from the prism side is indicated by a "O". Reflection and refraction over wavelength will differ when light input is applied from the opposite side of the prism.
 - ▶ Approximately 10% to 15% of absorption occurs in hybrid coating due to the properties in metallic coating reducing the transmitted or reflected light.
 - ▶ Phase retardation of light input will not be preserved. Use a wave-plate for phase compensation.
 - ▶ Wavelength dispersion of transmitted and reflected light is derived from refractive index and glass thickness and when diverging or introducing a focusing beam, chromatic aberration or spherical aberration may occur.

Specifications				
Part Number	Wavelength Range [nm]	A=B=C [mm]	Transmittance [%]	Polarization dependency Tp-Ts [%]
HBCH-10-550	400 - 700	10	45±10 (550nm)	<10
HBCH-15-550	400 - 700	15	45±10 (550nm)	<10
HBCH-20-550	400 - 700	20	45±10 (550nm)	<10
HBCH-10-NIR	700 - 1100	10	47±10 (900nm)	<20 (<10: 800 - 1100nm)
HBCH-15-NIR	700 - 1100	15	47±10 (900nm)	<20 (<10: 800 - 1100nm)
HBCH-20-NIR	700 - 1100	20	47±10 (900nm)	<20 (<10: 800 - 1100nm)
HBCH-10-IR	1300 - 1550	10	45±10 (1400nm)	<10
HBCH-15-IR	1300 - 1550	15	45±10 (1400nm)	<10
HBCH-20-IR	1300 - 1550	20	45±10 (1400nm)	<10

Typical Transmittance Data




Compatible Optic Mounts

PLH-25, -40 / KKD-25PHRO, -40PHRO


- Application Systems
- Optics & Optical Coatings
- Opto-Mechanics
- Bases
- Manual Stages
- Actuators & Adjusters
- MotORIZED Stages
- Light Sources & Laser Safety
- Index
- Guide
- Mirrors
- Beamsplitters
- Polarizers
- Lenses
- Multi-Element Optics
- Filters
- Prisms
- Substrates/Windows
- Optical Data
- Maintenance
- Selection Guide
- Half Mirror Cube
- Half Mirror Plate
- Application Note
- Beamsplitters
- Harmonic Separator
- Beam Samplers
- Others

Chromium cube half mirrors consist of two right angle prisms. One of them is coated with chromium (Cr) on the hypotenuse face. Half mirror divides input beam to reflectance and transmittance in 1:1. A beamsplitter of R:T=1:1 is called "Half Mirror".

- Four surfaces of the cube are coated with multi-layer anti-reflection coatings.
- Approximately one third of the input beam is lost because of absorption of chromium coating on the hypotenuse. However these beamsplitters are not wavelength, polarization and incident angle of the input beam dependent and therefore provide a highly neutral reflectivity.
- For cube beamsplitters, unlike plate beamsplitters, transmission beam deviations and ghosts rarely occur.




Schematic



Outline Drawing

(in mm)



- Tolerance
- A ±0.2
- B ±0.2
- C ±0.1

Specifications

Part Number	Wavelength Range [nm]	A=B=C [mm]
CSCH-10-550	400 – 700	10
CSCH-15-550	400 – 700	15
CSCH-20-550	400 – 700	20
CSCH-25-550	400 – 700	25
CSCH-30-550	400 – 700	30
CSCH-40-550	400 – 700	40
CSCH-50-550	400 – 700	50
CSCH-10-800	750 – 850	10
CSCH-15-800	750 – 850	15
CSCH-20-800	750 – 850	20

Specifications

Material	BK7
Surface flatness of substrate	$\lambda/4$
Beam Deviation	<5'
Coating	Hypotenuse surface: Chromium Four surfaces: Multi-layer anti-reflection coating
Incident angle	0°
Transmittance	Average 28±5% (The average value of the P-Polarization and the S-Polarization)
Divergence ratio (reflectance : transmittance)	1 : 1
Laser Damage Threshold	0.3J/cm ² (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	40-20
Clear aperture	85% of actual aperture

Guide

- ▶ Please contact our Sales Team for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division.


Attention

- ▶ Input beam from the prism side is indicated by a "O".
- ▶ Phase retardation of light input will not be preserved. Use a wave-plate for phase compensation.
- ▶ Wavelength dispersion of transmitted and reflected light is derived from refractive index and glass thickness and when diverging or introducing a focusing beam, chromatic aberration or spherical aberration may occur.
- ▶ The transmittance curves are based on actual measurements and may vary with manufacturing lots.
- ▶ The surface flatness is the reflected wavefront distortion of the surface before coating.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.


Typical Reflectance Data

T: Transmission

CSCH-550



CSCH-800



Compatible Optic Mounts


PLH-25, -40 / KKD-25PHRO, -40PHRO

Dielectric Cube Half Mirrors | CSMH


RoHS

Dielectric cube half mirrors consist of two right angle prisms. One of them is coated with dielectric multi-layer partial reflection coating on the hypotenuse face.


- Half mirror divides input beam into reflectance and transmittance at a 1:1 ratio. A beamsplitter with R:T (1:1 ratio) is called "Half Mirror".
- Four surfaces of the cube are coated with multi-layer anti-reflection coatings.
- The loss of input beam is minimized as there is no absorption from dielectric coating. However the reflection to transmittance ratio of these dielectric cube half mirrors vary depending on wavelength, polarization and the incident angle of input beam.



Schematic



Outline Drawing



- Tolerance
A ± 0.2
B ± 0.2
C ± 0.1

(in mm)

Specifications

Material	BK7
Surface flatness of substrate	$\lambda/4$
Beam Deviation	$<5'$
Coating	Hypotenuse surface: Dielectric multi-layer coating Four surfaces: Multi-layer anti-reflection coating
Incident angle	0°
Divergence ratio (reflectance : transmittance)	1 : 1
Polarization of the incident beam	Unpolarized light or 45 degrees Linear polarization or circular polarization
Laser Damage Threshold	$0.3J/cm^2$ (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	20-10
Clear aperture	85% of circle to actual dimension (80% of actual aperture for 5 and 7mm dimension (A=B=C) products.)

Guide

- ▶ Please contact our Sales Team for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division.

Attention

- ▶ Input beam from the prism side is indicated by a "O". Reflection and refraction over wavelength will differ when light input is applied from the opposite side of the prism.
- ▶ The transmittance curves are based on actual measurements and may vary with manufacturing lots.
- ▶ The surface flatness is the reflected wavefront distortion of the surface before coating.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.

Specifications

Part Number	Wavelength Range [nm]	A=B=C [mm]	Transmittance (The average value of the P-Polarization and the S-Polarization) [%]
CSMH-10-405	390 - 410	10	Average 50 ± 3
CSMH-12.7-405	390 - 410	12.7	Average 50 ± 3
CSMH-15-405	390 - 410	15	Average 50 ± 3
CSMH-20-405	390 - 410	20	Average 50 ± 3
CSMH-25-405	390 - 410	25	Average 50 ± 3
CSMH-30-405	390 - 410	30	Average 50 ± 3
CSMH-05-550	400 - 700	5	Average 50 ± 5
CSMH-07-550	400 - 700	7	Average 50 ± 5
CSMH-10-550	400 - 700	10	Average 50 ± 5
CSMH-12.7-550	400 - 700	12.7	Average 50 ± 5
CSMH-15-550	400 - 700	15	Average 50 ± 5
CSMH-20-550	400 - 700	20	Average 50 ± 5
CSMH-25-550	400 - 700	25	Average 50 ± 5
CSMH-30-550	400 - 700	30	Average 50 ± 5
CSMH-40-550	400 - 700	40	Average 50 ± 5
CSMH-50-550	400 - 700	50	Average 50 ± 5
CSMH-10-800	750 - 850	10	Average 50 ± 5
CSMH-12.7-800	750 - 850	12.7	Average 50 ± 5
CSMH-15-800	750 - 850	15	Average 50 ± 5
CSMH-20-800	750 - 850	20	Average 50 ± 5
CSMH-25-800	750 - 850	25	Average 50 ± 5
CSMH-30-800	750 - 850	30	Average 50 ± 5
CSMH-10-1400	1300 - 1550	10	Average 50 ± 5
CSMH-12.7-1400	1300 - 1550	12.7	Average 50 ± 5
CSMH-20-1400	1300 - 1550	20	Average 50 ± 5

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator


Beam Samplers

Others


Typical Transmittance Data

T: Transmission


CSMH-405




CSMH-550



CSMH-800



CSMH-1400



Compatible Optic Mounts

PLH-25, -40 / KKD-25PHRO, -40PHRO, -60PHRO / SHA-60RO

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others


We offer Half-Mirror optics designed for use in Ultraviolet, Visible and Infrared wavelengths. They can be used for both transmission and divergence of multi-wavelength lasers and white light sources.

Ultra broadband half-mirrors are used for spectrometry applications.


- PMH series have 4 types of ultra-broadband optics suitable for UV to IR applications.
- PSMH series have 3 types of ultra-broadband optics for ranges from Visible to NIR, which are used for optical communication applications.
- Dielectric multi-layer coated optics are an excellent choice for beam deviation applications because of its low absorption capabilities.
- Its low polarization characteristic can also be applied in beam deviation with a linear polarization laser or a laser light.
- Sigma Koki produces plate from optics that are light weight and maintains low dispersion with less aberration.
- Both wedge and plate type mirrors are made to have "low ghosting and low interference effects".



Schematic



Outline Drawing



Specifications

Material	BK7, Synthetic fused silica
Surface Flatness	$\lambda/10$
Coating	Front surface: Dielectric multi-layer coating Rear surface: Multi-layer anti-reflection coating
Incident angle	45°
Divergence ratio (reflectance : transmittance)	1 : 1
Surface Quality (Scratch-Dig)	10-5
Clear aperture	90% of actual aperture

Guide

- ▶ For customization, we can offer different sizes, wavelengths and deviation ratios. [Reference](#) B068
Please contact our Sales Division with your requests.
- ▶ For guaranteed higher reflectance accuracy and higher transmittance optics, please contact us.
- ▶ An arrow mark will be printed on the thick side of the wedge plate to indicate the surface of the mirror.

Attention

- ▶ When applying laser linear polarized light, the direction of polarization may affect the amount of reflectance and transmittance. For a divergence usage of 1:1 ratio, ensure the direction of polarization is set to 45 degrees or use a circular polarizer.
- ▶ When a laser light transmits through the optics, the light path may shift by a few millimeters horizontally due to refraction and the thickness of the wedge plate.
- ▶ The transmittance wavelength properties may be different if the incident angle is other than 45 degrees.
- ▶ Please check the arrow mark on the side of the wedge plate that indicates the coated surface.
- ▶ The phase difference of incident light cannot be preserved on light transmittance and reflectance; Please use a wave plate to compensate.

Ultra broadband

Part Number	Wavelength Range [nm]	Diameter ϕD [mm]	Thickness t [mm]	Material	Parallelism W	Transmittance (The average value of the P-Polarization and the S-Polarization) [%]	Laser Damage Threshold* [J/cm ²]
PMH-25.4C03-10-25/7	250 - 700	$\phi 25.4$	3	Synthetic fused silica	<5"	Average 50±10	0.5
PMH-30C03-10-25/7	250 - 700	$\phi 30$	3	Synthetic fused silica	<5"	Average 50±10	0.5
PMH-50C05-10-25/7	250 - 700	$\phi 50$	5	Synthetic fused silica	<5"	Average 50±10	0.5
PMH-25.4C03-10-3/10	300 - 1000	$\phi 25.4$	3	Synthetic fused silica	<5"	Average 50±10	0.5
PMH-30C03-10-3/10	300 - 1000	$\phi 30$	3	Synthetic fused silica	<5"	Average 50±10	0.5
PMH-50C05-10-3/10	300 - 1000	$\phi 50$	5	Synthetic fused silica	<5"	Average 50±10	0.5
PMH-25.4C03-10-6/18	600 - 1800	$\phi 25.4$	3	BK7	<5"	Average 50±10	0.5
PMH-30C03-10-6/18	600 - 1800	$\phi 30$	3	BK7	<5"	Average 50±10	0.5
PMH-50C05-10-6/18	600 - 1800	$\phi 50$	5	BK7	<5"	Average 50±10	0.5
PMH-25.4C03-10-4/20	400 - 2000	$\phi 25.4$	3	BK7	<5"	Average 50±10	0.5
PMH-30C03-10-4/20	400 - 2000	$\phi 30$	3	BK7	<5"	Average 50±10	0.5
PMH-50C05-10-4/20	400 - 2000	$\phi 50$	5	BK7	<5"	Average 50±10	0.5

* Laser pulse width 10ns, repetition frequency 20Hz

Compatible Optic Mounts

BHAN-30S, -50S / MHG-HS25-NL, MP30-NL, MP50-NL

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator


Beam Samplers

Others


Typical Transmittance Data

T: Transmission


PMH-25/7




PMH-3/10



PMH-6/18



PMH-4/20



Broadband


Part Number	Wavelength Range [nm]	Diameter ϕ D [mm]	Thickness t [mm]	Material	Parallelism W	Transmittance (The average value of the P-Polarization and the S-Polarization) [%]	Laser Damage Threshold* [J/cm ²]
PSMH-25.4C03-10-550	400 – 700	ϕ 25.4	3	BK7	<5"	Average 50 \pm 5	2.1
PSMH-30C03-10-550	400 – 700	ϕ 30	3	BK7	<5"	Average 50 \pm 5	2.1
PSMH-30C05-10W-550	400 – 700	ϕ 30	5	BK7	1° \pm 5'	Average 50 \pm 5	2.1
PSMH-40C04-10-550	400 – 700	ϕ 40	4	BK7	<5"	Average 50 \pm 5	2.1
PSMH-50C05-10-550	400 – 700	ϕ 50	5	BK7	<5"	Average 50 \pm 5	2.1
PSMH-50C08-10W-550	400 – 700	ϕ 50	8	BK7	1° \pm 5'	Average 50 \pm 5	2.1
PSMH-30C03-10-800	700 – 900	ϕ 30	3	BK7	<5"	50 \pm 3 (800nm)	2.1
PSMH-30C05-10W-800	700 – 900	ϕ 30	5	BK7	1° \pm 5'	50 \pm 3 (800nm)	2.1
PSMH-50C05-10-800	700 – 900	ϕ 50	5	BK7	<5"	50 \pm 3 (800nm)	2.1
PSMH-50C08-10W-800	700 – 900	ϕ 50	8	BK7	1° \pm 5'	50 \pm 3 (800nm)	2.1
PSMH-30C03-10-1400	1300 – 1550	ϕ 30	3	BK7	<5"	50 \pm 3 (1400nm)	2.1
PSMH-30C05-10W-1400	1300 – 1550	ϕ 30	5	BK7	1° \pm 5'	50 \pm 3 (1400nm)	2.1

* Laser pulse width 10ns, repetition frequency 20Hz


Typical Transmittance Data

T: Transmission


PSMH-550



PSMH-800



PSMH-1400



Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others

Thin Plate Beamsplitter | MPSMH

RoHS



Catalog Code

W3017


Extremely thin beamsplitter.

It can be inserted into an optical light path without any beam shift or chromatic dispersion for any light transmittance application.

- We offer 2 choices of thickness, 300um and 90um.
- Dielectric multi-layer optical coating with reflectance and transmittance ratios at 1:1.
- Dielectric multi-layer optical coating on the front surface and AR coating on the rear surface to provide a mirror with no loss of power.
- Due to the fabrication method, these offer good durability and high resistance against vibration making them an excellent alternative to traditional pellicle beamsplitters.

**Schematic****Outline Drawing**

(in mm)

**Specifications**

Part Number	Wavelength Range [nm]	Optics Thickness [mm]	Surface flatness after coating
MPSMH-30C0.3-1-550	400 – 700	0.3±0.03	Reflectance: λ Transmittance: λ
MPSMH-30C0.09-1-550	400 – 700	0.09±0.01	Reflectance: Polishing Transmittance: Polishing

Specifications

Material	Synthetic fused silica
Coating	Front surface: Dielectric multi-layer coating Rear surface (45 degrees taper hole): Anti-reflection coating
Incident angle	45°
Transmittance	Average 50±5% (The average value of the P-Polarization and the S-Polarization)
Divergence ratio (reflectance : transmittance)	1 : 1
Surface Quality (Scratch-Dig)	40-20
Clear aperture	ϕ 10mm
Frame specification	Frame glass: Synthetic fused silica Outer frame: Aluminum Finishing: Matt black almite

Guide


- ▶ For customization, we can offer different sizes, wavelengths and deviation ratios. [Reference](#) B068
Please contact our Sales Division with your requests.

Attention

- ▶ Thin beamsplitters are extremely thin and fragile. Special care must be taken during cleaning and handling.
- ▶ When removing dust from the surface, do not use optics tissue paper to clean. Use a compress gas spray instead.
- ▶ When applying a laser linear polarized light, the direction of polarization may affect the amount of reflectance and transmittance. For divergence usage of 1:1 ratio, ensure the direction of polarization is set to 45 degrees or use a circular polarizer.
- ▶ The transmittance wavelength properties may be different if the incident angle is other than 45 degrees.
- ▶ Avoid pushing the glass retainer as the mirror can bend or break. When handling, please use the other metal frame.
- ▶ The surface reflectance accuracy may deteriorate when used outside recommended operating temperature.
- ▶ The phase difference of incident light cannot be preserved on light transmittance and reflectance. Please use a wave plate to compensate.

Typical Transmittance Data

T: Transmission

**Compatible Optic Mounts**

MHG-HS30-NL / BHAN-30S

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters


Harmonic Separator

Beam Samplers


Others

Laser line plate mirrors are plate beamsplitters that are optically coated with dielectric multi-layer on the front surface of optical parallels or wedged substrates. The rear surface is coated with multi-layer anti-reflection.


- Half mirrors divide input beam into reflectance and transmittance ratio of 1:1. A beamsplitter of R:T=1:1 is called "Half Mirror".
- Any loss from the input beams of this product is minimized because dielectric coatings have no absorption properties. However, the input ratio of reflection to transmission depends on wavelength, polarization and angle of incident of input beam.
- Plate beamsplitters have beam deviations on transmission and ghost on rear surface reflections. Wedged substrates are used to prevent ghosting.




Schematic



Optical Parallel



Wedged Substrate



Specifications

Material	BK7, Synthetic fused silica, CaF ₂
Surface Flatness	$\lambda/10$ (PSMH-157 is Polished)
Coating	Front surface: Dielectric multi-layer partial reflection coating Rear surface: Multi-layer anti-reflection coating
Incident angle	45°
Divergence ratio (reflectance : transmittance)	1 : 1
Surface Quality (Scratch-Dig)	10-5 (PSMH-157: 40-20)
Clear aperture	90% of actual aperture

Guide


- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ We also have ultra-wideband, broadband and cube types.
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division with your requests.
- ▶ Wedged types are marked with an arrow on the side of the substrate indicating the thickest point of the wedge.

Attention

- ▶ Should these products not function as a half mirror, please check the polarization characteristics of the light source. Do note that LD laser is linear in polarization.
- ▶ The beam deviation at transmission of a wedged beamsplitter is large compared to a one made of optical parallel.
- ▶ The amount of beam deviation of a beamsplitter depends on the thickness of the substrate and the wavelength or the incident angle of the input beam.
- ▶ Transmission curves are based on actual measurements and may vary with manufacturing lots.
- ▶ Surface flatness is the reflected wavefront distortion of the surface prior to coating.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.

Outline Drawing

(in mm)



●Tolerance
Diameter $\phi D_{\pm 0.1}$
Thickness $t \pm 0.1$

Compatible Optic Mounts

BHAN-30S, -50S / MHG-MP30-NL, MP50-NL

Laser Line							
Part Number	Wavelength Range [nm]	Diameter ϕ D [mm]	Thickness t [mm]	Material	Parallelism W	Reflectance:R Transmittance:T (The average value of the P-Polarization and the S-Polarization) [%]	Laser Damage Threshold* [J/cm ²]
PSMH-30C03-P-157	157	ϕ 30	3	CaF ₂	<3'	R=40±10	0.5
PSMH-50C05-P-157	157	ϕ 50	5	CaF ₂	<3'	R=40±10	0.5
PSMH-30C03-10-193	193	ϕ 30	3	Synthetic fused silica	<5"	T=45±5	1
PSMH-30C05-10W-193	193	ϕ 30	5	Synthetic fused silica	1°±5'	T=45±5	1
PSMH-50C05-10-193	193	ϕ 50	5	Synthetic fused silica	<5"	T=45±5	1
PSMH-50C08-10W-193	193	ϕ 50	8	Synthetic fused silica	1°±5'	T=45±5	1
PSMH-30C03-10-248/266	248 – 266	ϕ 30	3	Synthetic fused silica	<5"	T=50±3	2
PSMH-30C05-10W-248/266	248 – 266	ϕ 30	5	Synthetic fused silica	1°±5'	T=50±3	2
PSMH-50C05-10-248/266	248 – 266	ϕ 50	5	Synthetic fused silica	<5"	T=50±3	2
PSMH-50C08-10W-248/266	248 – 266	ϕ 50	8	Synthetic fused silica	1°±5'	T=50±3	2
PSMH-30C03-10-308/355	308 – 355	ϕ 30	3	Synthetic fused silica	<5"	T= Average 50±5	2
PSMH-30C05-10W-308/355	308 – 355	ϕ 30	5	Synthetic fused silica	1°±5'	T= Average 50±5	2
PSMH-50C05-10-308/355	308 – 355	ϕ 50	5	Synthetic fused silica	<5"	T= Average 50±5	2
PSMH-50C08-10W-308/355	308 – 355	ϕ 50	8	Synthetic fused silica	1°±5'	T= Average 50±5	2
PSMH-30C03-10-405	390 – 410	ϕ 30	3	BK7	<5"	T=50±3	2.1
PSMH-30C05-10W-405	390 – 410	ϕ 30	5	BK7	1°±5'	T=50±3	2.1
PSMH-50C05-10-405	390 – 410	ϕ 50	5	BK7	<5"	T=50±3	2.1
PSMH-50C08-10W-405	390 – 410	ϕ 50	8	BK7	1°±5'	T=50±3	2.1
PSMH-30C03-10-1064	1064	ϕ 30	3	BK7	<5"	T=50±3	20
PSMH-30C05-10W-1064	1064	ϕ 30	5	BK7	1°±5'	T=50±3	20
PSMH-50C05-10-1064	1064	ϕ 50	5	BK7	<5"	T=50±3	20
PSMH-50C08-10W-1064	1064	ϕ 50	8	BK7	1°±5'	T=50±3	20

*Laser pulse width 10ns (PSMH-157: 20ns), repetition frequency 20Hz

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers


Others




Typical Reflectance Data & Typical Transmittance Data

R: Reflectance T: Transmission


PSMH-157




PSMH-193




PSMH-248/266




PSMH-308/355



PSMH-405



PSMH-1064



Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others


Chromium plate half mirrors are plate beamsplitters that are coated with chromium (Cr) on the front surface of optical parallels or wedged substrates.

The other surface is coated with multi-layer anti-reflection.

- Half mirror divides input beam into reflectance and transmittance ratio of 1:1. A beamsplitter of R:T=1:1 is called "Half Mirror".
- Approximately one third of the input beam is lost because of the absorption of chromium. However these beamsplitters do not depend on wavelength, polarization and angle of incidence of the input beam, and provide a highly neutral reflectivity.
- Plate beamsplitters have beam deviations on transmission and ghost on rear surface reflections. Wedged substrates are used to prevent ghosting.




Schematic



Outline Drawing

(in mm)



Specifications

Part Number	Wavelength Range [nm]	Diameter ϕD [mm]	Thickness t [mm]	Parallelism W
PSCH-25.4C03-10-550	400 - 700	$\phi 25.4$	3	$<5''$
PSCH-30C03-10-550	400 - 700	$\phi 30$	3	$<5''$
PSCH-30C05-10W-550	400 - 700	$\phi 30$	5	$1'' \pm 5''$
PSCH-40C04-10-550	400 - 700	$\phi 40$	4	$<5''$
PSCH-50C05-10-550	400 - 700	$\phi 50$	5	$<5''$
PSCH-50C08-10W-550	400 - 700	$\phi 50$	8	$1'' \pm 5''$
PSCH-60C06-10-550	400 - 700	$\phi 60$	6	$<5''$
PSCH-100C10-10-550	400 - 700	$\phi 100$	10	$<5''$
PSCH-100C15-10W-550	400 - 700	$\phi 100$	15	$1'' \pm 5''$
PSCH-25.4C03-10-800	750 - 850	$\phi 25.4$	3	$<5''$
PSCH-30C03-10-800	750 - 850	$\phi 30$	3	$<5''$
PSCH-30C05-10W-800	750 - 850	$\phi 30$	5	$1'' \pm 5''$
PSCH-50C05-10-800	750 - 850	$\phi 50$	5	$<5''$
PSCH-50C08-10W-800	750 - 850	$\phi 50$	8	$1'' \pm 5''$

Compatible Optic Mounts

BHAN-30S, -50S / MHAN-25.4S, -40S, -60S / MHG-MP25-NL, MP30-NL, MP50-NL

Specifications

Material	BK7
Surface Flatness	$\lambda/10$
Coating	Front surface: Chromium Rear surface: Multi-layer anti-reflection coating
Incident angle	45°
Transmittance	Average $30 \pm 5\%$ (The average value of the P-Polarization and the S-Polarization)
Divergence ratio (reflectance : transmittance)	1 : 1
Laser Damage Threshold	0.25 J/cm^2 (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	40-20
Clear aperture	90% of actual aperture


Guide

- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division with your requests.


Attention

- ▶ The beam deviation at transmission of a wedged beamsplitter is large compared to a one made of optical parallel.
- ▶ The amount of beam deviation of a beamsplitter depends on the thickness of the substrate, the wavelength or the angle of incidence of the input beam.
- ▶ Transmission curves are based on actual measurements and may vary with manufacturing lots.
- ▶ Surface flatness is the reflected wavefront distortion of the surface prior to coating.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.

Optical Parallel




Wedged Substrate




Typical Transmittance Data

T: Transmission

PSCH-550



PSCH-800



About light behaviour on a beamsplitter

A half mirror is designed with reflectance and transmission of light with a 1:1 ratio. If light incident direction and polarization conditions change, it may impact the ratio.

Reflectance and transmittance properties of the incident light direction

Chrome coating and multi-wavelength coating application.

Reflection properties change when light is projected onto the coated and black surfaces. Any configuration similar to Michelson interferometer may require both sides to have incident light. In this case, light ratios may be unbalanced.



Choose the following set up if the light incident direction can be selected.

Incident light onto the coated surface of plate type beamsplitter.


Incident light onto the \odot mark surface for cube type beamsplitter.

If the Incident light is on the wrong surface, the specifications mentioned in the catalogue cannot be realized.

Comparison reflectance and transmittance properties of the incident light direction in the chromium plate half mirror.



The difference in reflectance due to the incident direction occurs when there is absorption in the coating. It does not occur in the dielectric multilayer coating.




The reflectance and the transmittance of a polarized light incident

In case fo using Laser

Light emitted from the laser is linearly polarized light. Because of this, even though it is used in the experiments and the optical system which are not related to the polarization, it is necessary to take into account the polarization characteristics of the beam splitter.

The transmittance and the reflectance may change in accordance with the type of beamsplitter and its polarization direction.

To split the light into a balanced light ratio, a non-polarized beam splitter (NPCH) is recommended. The polarization properties of the laser has no influence to it.



Dielectric Cube Beamsplitters

CSM


RoHS

Catalog Code W3020


Cube beamsplitters with dielectric multi-layer coated on the hypotenuse face of a 45° right angle prism.

Divides beams at reflected (R): transmission light (T) ratio of 1:2 or 1:3.


- Anti-reflection coating (AR coat) is applied to the incident and outgoing faces.
- The dielectric multi-layer films has virtually zero light absorption and very low light intensity loss. However, transmittance and reflectance may change according to wavelength, polarization and incident angles.
- In contrast to plate type half mirrors, cube mirrors have no ghosting or transmission optical path deviation.



Schematic



Outline Drawing (in mm)



Specifications

Material	BK7
Surface Flatness	$\lambda/4$
Wavelength Range	400 – 700nm
Beam Deviation	<5'
Coating	Hypotenuse surface: Dielectric multi-layer coating Four surfaces: Multi-layer anti-reflection coating
Incident angle	0°
Polarization of the incident beam	Unpolarized light or 45 degrees Linear polarization or circular polarization
Laser Damage Threshold	0.3J/cm ² (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	20-10
Clear aperture	85% of actual aperture

Guide

- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division with your requests.



Attention

- ▶ Input beam from the prism side is indicated by a "○" (hypotenuse coated side).
- ▶ The transmission curve on the graph is based on actual measurements and may vary from different production lots.
- ▶ Phase retardation of inputting light will not be preserved. Please use waveplate for phase compensation.
- ▶ Use only non-polarized light or circular polarized light as incident light for dielectric multi-layer coated beam splitters. Using polarized light may result in R:T ratios that vary according to polarization components.
- ▶ Dielectric multi-layer coated cube half mirrors sometimes do not function effectively. If this should occur, first check the polarization characteristics of the light source (laser) and keep in mind that lasers used in the semiconductor field emit a linear polarized light.

Specifications

Part Number	Reflectance : Transmittance	A=B=C [mm]	Transmittance at 550nm	Transmittance at 400-700nm
			(The average value of the P-Polarization and the S-Polarization) [%]	(The average value of the P-Polarization and the S-Polarization) [%]
CSM33-10-550	1 : 2	10	67±5	<80
CSM33-20-550	1 : 2	20	67±5	<80
CSM25-10-550	1 : 3	10	75±5	<90
CSM25-20-550	1 : 3	20	75±5	<90

Typical Transmittance Data




Compatible Optic Mounts


PLH-25, -40 / KKD-25PHRO, -40PHRO

Plate-type beamsplitters are dielectric multi-layer coatings on a parallel plate or a wedge substrate. They are designed to divide beams at a reflected light (R) : transmission light (T) ratio of 1:2 or 1:3. The rear surface is coated with anti-reflection (AR).

- The dielectric multi-layer films have virtually zero light absorption and very low light intensity loss. However, transmittance and reflectance may vary according to wavelength, polarization and incident angles. Some deviation of the transmission optical path or ghosting may occur. To prevent ghosting, use wedge beamsplitters.




Schematic



Outline Drawing

(in mm)



Specifications

Material	BK7
Surface Flatness	$\lambda/10$
Coating	Front surface: Dielectric multi-layer coating Rear surface: Multi-layer anti-reflection coating
Wavelength Range	400 – 700nm
Incident angle	45°
Laser Damage Threshold	2.1J/cm ² (Laser pulse width 10ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	10-5
Clear aperture	90% of actual aperture

Guide

- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division with your requests.
- ▶ Wedged types are marked with an arrow on the side of the substrate indicating the thickest point of the wedge.

Attention

- ▶ The transmission curve on the graph is based on actual measurements and may vary from different production lots.
- ▶ Surface flatness is the reflected wavefront distortion of the surface prior to coating.
- ▶ Compared to precision parallel plate type splitters, wedged substrate type beam splitters can prevent ghosting caused by rear surface reflection and significantly increase the displacement of the optical path.
- ▶ Dielectric multi-layer coated beamsplitters sometimes do not function effectively in specified R:T ratios. If this should occur, first check the polarization characteristics of the light source (laser). Do keep in mind that lasers used in the semiconductor field emit a linear polarized light.
- ▶ Use only non-polarized light or circular polarized light as incident light for dielectric multi-layer coated beam splitters. Using polarized light may result in R:T ratios that vary according to polarization components.


Specifications

Part Number	Reflectance : Transmittance	Diameter φD [mm]	Thickness t [mm]	Parallelism W	Transmittance at 550nm (The average value of the P-Polarization and the S-Polarization)		Transmittance at 400-700nm (The average value of the P-Polarization and the S-Polarization)	
					[%]	[%]	[%]	[%]
PSM33-25.4C03-10-550	1 : 2	φ25.4	3	<5"	67±3	<80		
PSM33-30C03-10-550	1 : 2	φ30	3	<5"	67±3	<80		
PSM33-30C05-10W-550	1 : 2	φ30	5	1°±5'	67±3	<80		
PSM25-25.4C05-10-550	1 : 3	φ25.4	3	<5"	75±3	<90		
PSM25-30C03-10-550	1 : 3	φ30	3	<5"	75±3	<90		
PSM25-30C05-10W-550	1 : 3	φ30	5	1°±5'	75±3	<90		


Typical Transmittance Data

T: Transmission

PSM33



PSM25



Compatible Optic Mounts

BHAN-30S / MHAN-25.4DS / MHG-MP25-NL, MP30-NL

Variable Beamsplitter Light path corrector


VBS
WSQNA/WBNA

RoHS
RoHS


With a variable beamsplitter, the incident angle of a laser can be changed. The (R:T) ratios can also be modified.

This is commonly used to adjust the light intensity of the laser without a variable adjustment of the light intensity or the laser to be stabilized.

- The variable beamsplitter has a dielectric multilayer coating which has excellent durability and light resistance.
- The beam shift caused by the tilt of the beamsplitter can be removed by using a correcting plate. (See how to use)
- It can be used for arbitrary polarization. However, the transmittance characteristic depends on the polarization.




Schematic



Outline Drawing

(in mm)



Specifications

Material	BK7, Synthetic fused silica
Surface Flatness	λ
Parallelism	<5"
Coating	VBS Front surface: Dielectric multi-layer Coating Rear surface: Multi-layer anti-reflection coating WBMA, WSQMA Both surfaces: Multi-layer anti-reflection coating
Surface Quality (Scratch-Dig)	10-5
Clear aperture	Circle that internally connected to 90% of the side length
Effective beam incident diameter	Ellipsoidal 30x43mm (Angle of inclination)

Guide

▶ Different size, wavelength and deviation not mentioned on-line or in our catalog are available as custom product upon request.

Reference B068

▶ We offer the most comprehensive range of beamsplitter holders and stages to select from. Let us know the angle of your choice.

▶ This variable attenuator (model SVAB) can be used in a system and is available.



Attention

▶ When using with high power laser, make sure to execute at the end edge of the reflected light.

▶ The reflectance properties of the optics may change in a high temperature environment.

▶ When adjusting the transmittance, the incident angle may change and cause the light path to shift. To correct this, please use the light path corrector (model WSQNA/WBNA)

▶ For a large beam size of 30mm diameter or larger and used it at a high inclination level, the beam can be cut at the reflected area.

▶ For "P" polarization use, make sure that the incident angle is at 45 degrees or more.

Variable beamsplitter

Part Number	Wavelength Range [nm]	Transmittance of S polarization ($\theta=0^\circ$) [%]	Transmittance of S polarization ($\theta=45^\circ$) [%]	Material	Laser Damage Threshold* [J/cm ²]
VBS-50S03-1-266	266	>90	<5	Synthetic fused silica	1
VBS-50S03-1-355	355	>93	<5	Synthetic fused silica	1
VBS-50S03-1-532	532	>95	<5	BK7	2.5
VBS-50S03-1-1064	1064	>95	<5	BK7	3.5

* Laser pulse width 10ns, repetition frequency 20Hz

Light path corrector

Part Number	Wavelength Range [nm]	Transmittance of S polarization ($\theta=0^\circ - 45^\circ$) [%]	Material	Laser Damage Threshold* [J/cm ²]
WSQNA-50S03-1-266-0/45D	266	Average 97	Synthetic fused silica	1
WSQNA-50S03-1-355-0/45D	355	Average 97	Synthetic fused silica	1
WBNA-50S03-1-532-0/45D	532	Average 98	BK7	2.5
WBNA-50S03-1-1064-0/45D	1064	Average 98	BK7	3.5

* Laser pulse width 10ns, repetition frequency 20Hz

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator


Beam Samplers

Others


Typical Transmittance Data

T: Transmission (S polarization)


VBS-266 / WSQNA-266




VBS-355 / WSQNA-355



VBS-532 / WBNA-532



VBS-1064 / WBNA-1064




Sample of use

The variable beam splitter can be used individually. When modifying the incident angle, optics thickness and its refractive properties, a shift may occur in the light path. To reduce this shift, we highly recommend a light path corrector. Please see image below.

- Place the variable beamsplitter onto a rotation stage to allow an angle adjustment.
- Install the light path corrector onto a rotating stage.
- Position the light path corrector at a similar angle with the variable beamsplitter on an opposite side.
- If the reflected light of the variable beamsplitter is not used, make sure to place a light cut-off material or a beam diffuser at the edge-end of the light.
- The power of the reflected light from the light path corrector must be cut off at the edge-end of the light.

For part structure, please contact our International Sales Division.




Compatible Optic Mounts


CHA-60, -60F

Harmonic separators are designed to separate specific YAG harmonics from other harmonics. We offer three standard wavelength (1064nm, 532nm, 355nm) reflectance YAG harmonics.

- These mirrors are coated with multi-layered dielectric with different refractive index, using BK7 optical parallels with $\lambda/10$ surface flatness and parallelism of 5 arc second. The other surface is coated with multi-layer anti-reflection.
- These mirrors are used at 45° incident angle to reflect specific wavelength beam and transmits other YAG wavelengths.
- For plate type, you can use a large laser beam diameter.




Schematic



Outline Drawing

(in mm)



Specifications

Material	BK7
Surface Flatness	$\lambda/10$
Coating	Front surface: Dielectric multi-layer coating Rear surface: Multi-layer anti-reflection coating
Angle of Incidence	45°
Parallelism	<5"
Surface Quality (Scratch-Dig)	10-5
Clear aperture	90% of actual aperture

Guide

- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) [Reference](#) B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Team with your requests.

Attention

- ▶ The reflection surface is indicated with an arrow on the side of substrate.
- ▶ The reflectance curves are based on actual measurements and may vary from different manufacturing lots.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.
- ▶ The reflectance in the specifications list is at random polarization or (p-polarization reflectance + s-polarization reflectance) / 2.

For Reflected wavelength : 355nm, Transmitted wavelength : 532, 1064nm

Part Number	Diameter ϕD [mm]	Thickness t [mm]	Reflectance at 355nm (The average value of the P-Polarization and the S-Polarization)		Transmittance at 532·1064nm (The average value of the P-Polarization and the S-Polarization)		Laser Damage Threshold* [J/cm ²]
			[%]	[%]	[%]	[%]	
YHS-25.4C05-355	$\phi 25.4$	5	>99.5	>99.5	>85	>85	5
YHS-30C05-355	$\phi 30$	5	>99.5	>99.5	>85	>85	5
YHS-50C08-355	$\phi 50$	8	>99.5	>99.5	>85	>85	5

*Laser pulse width 10ns, repetition frequency 20Hz

For Reflected wavelength : 532nm, Transmitted wavelength : 1064nm

Part Number	Diameter ϕD [mm]	Thickness t [mm]	Reflectance at 532nm (The average value of the P-Polarization and the S-Polarization)		Transmittance at 1064nm (The average value of the P-Polarization and the S-Polarization)		Laser Damage Threshold* [J/cm ²]
			[%]	[%]	[%]	[%]	
YHS-25.4C05-532	$\phi 25.4$	5	>99.5	>99.5	>95	>95	8
YHS-30C05-532	$\phi 30$	5	>99.5	>99.5	>95	>95	8
YHS-50C08-532	$\phi 50$	8	>99.5	>99.5	>95	>95	8

*Laser pulse width 10ns, repetition frequency 20Hz

For Reflected wavelength : 1064nm, Transmitted wavelength : 532nm

Part Number	Diameter ϕD [mm]	Thickness t [mm]	Reflectance at 1064nm (The average value of the P-Polarization and the S-Polarization)		Transmittance at 532nm (The average value of the P-Polarization and the S-Polarization)		Laser Damage Threshold* [J/cm ²]
			[%]	[%]	[%]	[%]	
YHS-25.4C05-1064	$\phi 25.4$	5	>99.5	>99.5	>90	>90	20
YHS-30C05-1064	$\phi 30$	5	>99.5	>99.5	>90	>90	20
YHS-50C08-1064	$\phi 50$	8	>99.5	>99.5	>90	>90	20

*Laser pulse width 10ns, repetition frequency 20Hz

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator


Beam Samplers

Others


Typical Reflectance Data

R: Reflectance


YHS-355



YHS-532



YHS-1064



Compatible Optic Mounts

MHG-HS25-NL, HS30-NL / MHG-PM50-NL / BHAN-30S, -50S

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others

Beam Samplers | BS4

RoHS


Catalog Code W3025

A beam sampler behaves like a plate beamsplitter, it has the ability to reflect approximately 5.2% of the total beam.

- Uncoated surface of optical parallels or wedged substrates are used as reflection surfaces. The rear surfaces are coated with multi-layer anti-reflection.
- These products have transmitted beam deviation and ghosting of the rear surface reflections due to the characteristics of plate beamsplitters.
- Wedged beam samplers with AR coating on the rear surface should be selected to prevent ghosting.




Schematic



Outline Drawing


(in mm)



φ30·φ50

Part Number	Wavelength Range [nm]	Diameter ϕD [mm]	Thickness t [mm]	Parallelism W
BS4-25.4C03-10-550	400 - 700	$\phi 25.4$	3	$<5''$
BS4-30C03-10-550	400 - 700	$\phi 30$	3	$<5''$
BS4-30C05-10W-550	400 - 700	$\phi 30$	5	$1'' \pm 5'$
BS4-50C05-10-550	400 - 700	$\phi 50$	5	$<5''$
BS4-50C08-10W-550	400 - 700	$\phi 50$	8	$1'' \pm 5'$

Typical Reflectance Data (Calculated value)



Compatible Optic Mounts

BHAN-30S, -50S / MHG-MP25-NL, MP30-NL

Specifications

Material	BK7
Surface Flatness	$\lambda/10$
Coating	Front Surface: Uncoated Rear Surface: Visible multi-layer anti-reflection coating
Incident angle	45°
Divergence ratio (reflectance : transmittance)	5 : 95 (The average value of the P-Polarization and the S-Polarization)
Laser Damage Threshold	$4J/cm^2$ (Laser pulse width 4ns, repetition frequency 20Hz)
Surface Quality (Scratch-Dig)	10-5
Clear aperture	90% of actual aperture


Guide

- ▶ Please contact our Sales Division for customized products. (Customized on size, wavelength or R:T, etc.) **Reference** B068
- ▶ For a guarantee in reflected wavefront error or transmitted wavefront error, please contact our Sales Division with your requests.
- ▶ Wedged types are marked with an arrow on the side of the substrate indicating the thickest point of the wedge.


Attention

- ▶ The Beam Sampler has a reflectance of 5.2% when the material is BK7 and the input beam is unpolarized or circularly polarized.
- ▶ The transmitted beam deviation of a wedged beamsplitter is larger than with a beamsplitter made of an optical parallel.
- ▶ The amount of beam deviation of a beamsplitter depends on thickness of the substrate and the wavelength of the incident angle of the input beam.
- ▶ Be sure to wear laser safety goggles when checking optical path and adjusting optical axis.


Optical Parallel



Wedged Substrate




Reflectivity change in accordance with incident angle (Calculated value)




The polka dot beamsplitter is a beamsplitter that has aluminum coating of halftone dots (polka dots) on the glass substrate. It has a low dependence on the incident angle and can be used in a wide range of wavelengths from ultraviolet region to the infrared region.


- Reflectance to transmittance ratio is manufactured by adjusting the area ratio of the points that have been coated.
- Unlike dielectric beamsplitters, the polka dot beamsplitter reflectance and transmittance ratio does not change as the incident angle changes.
- There are two sizes available $\phi 25.4\text{mm}$ and $\phi 50.8\text{mm}$ diameter and three types of reflectance to transmittance ratio, 7:3, 5:5 and 3:7.



Schematic



Outline Drawing



Specifications

Material	Synthetic fused silica
Parallelism	<3'
Coating	Front Surface: Al+MgF ₂ Rear Surface: Uncoated
Recommended angle of incidence	0 – 45°
Wavelength range	250 – 2200nm
Surface Quality (Scratch-Dig)	80-50
Dot pitch	0.3mm
Clear aperture	Circle except surrounding 1.5mm


Guide

▶ We also offer different sizes, wavelengths and ratios that are not listed on our website or in our catalog. [Reference](#) B068

Attention

- ▶ When used with a laser beam with high interference, diffraction occurs.
 - ▶ When light is incident, scattering light by the halftone dot occurs.
 - ▶ By the effect of the refractive index and the thickness of the substrate, the optical path of the transmitted light over the incident light will move by 0.5 extent parallel.
 - ▶ When the incident beam diameter is very thin, it is not possible to separate into the split ratio.
 - ▶ Do not clean with water or solvents. It may cause surface deterioration.
 - ▶ Please use in the environments which are non-condensing and less dust.
- If the dust or dirt is deposited, please do not blow but blow it off gently with dried air.

Typical Transmittance Data



Specifications

Part Number	Reflectance : Transmittance	Diameter ϕD [mm]	Thickness t [mm]	Transmission (Wavelength Range 555nm, Angle of Incidence : 0°) [%]
PDBS70-25.4C1.5	70 : 30	$\phi 25.4$	1.5	30%
PDBS70-50.8C1.5	70 : 30	$\phi 50.8$	1.5	30%
PDBSH-25.4C1.5	50 : 50	$\phi 25.4$	1.5	50%
PDBSH-50.8C1.5	50 : 50	$\phi 50.8$	1.5	50%
PDBS30-25.4C1.5	30 : 70	$\phi 25.4$	1.5	70%
PDBS30-50.8C1.5	30 : 70	$\phi 50.8$	1.5	70%

Compatible Optic Mounts

P25-NL, MP50.8-NL / MHAN-25.4S, -50.8S

Contact sheet for Special Half Mirror and Beamsplitter and Beam Sampler Estimation Order

Date

To: Sigma Koki Co., Ltd. **FAX +81-3-5638-6550**

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate





Application Note


Beamsplitters

Harmonic Separator

Beam Samplers

Others

Affiliation (Organization Name)								
Department			Name					
TEL		FAX		E-mail				
Country Adress								
Name & Designation		(Tentative name is okay)						
Drawing Number			Estimate		<input type="checkbox"/> Yes: by Date		<input type="checkbox"/> No	
Desired Delivery Date			Budget		JP Yen			
Quantity		pieces						
Substrates <small>If you do not specify a dimension tolerance is outside the standard tolerance.</small>	Standard product						If you are using a substrate of standard product, please fill in the product number.	
	Custom-made	Material		<input type="checkbox"/> BK7 <input type="checkbox"/> Synthetic fused silica <input type="checkbox"/> Other ()				
			ϕA	mm	c	mm		
			a	mm	d	mm		
			b	mm	e	mm		
			t	mm				
	Surface flatness of substrate		(at $\lambda = 632.8\text{nm}$)		In the case of plate beam splitter, please specify surface flatness and parallelism and wedge of substrate.			
Parallelism		Wedge		°				
Wavelength Range		$\lambda =$ nm		Incident angle		$\theta =$ °		
Type of Coating	Metallic Coating		<input type="checkbox"/> Half Mirror <input type="checkbox"/> Beamsplitter		Type of Light Source			
			<input type="checkbox"/> Other ()					
	R : T		:		Power or Energy			
					pulse width			
					Repetition frequency			
Dielectric multi-layer coating		<input type="checkbox"/> Half Mirror <input type="checkbox"/> Beamsplitter		Beam Size			mm	
		<input type="checkbox"/> Dichroic Mirror						
		<input type="checkbox"/> Other ()						
R % T %		%		* There was a more detailed specification, please fill in this field.				
AR coat		<input type="checkbox"/> Multilayer antireflection coating (SLAR)						
		<input type="checkbox"/> Multilayer antireflection coating (MLAR)						
		<input type="checkbox"/> Other ()						
Polarization of the incident beam		<input type="checkbox"/> UnPolarization						
		<input type="checkbox"/> Circularly polarization or 45 degrees direction of lineraly polarization						
		<input type="checkbox"/> P-polarization <input type="checkbox"/> S-polarization						



Application
Systems

**Optics &
Optical
Coatings**

Opto-
Mechanics

Bases

Manual
Stages

Actuators &
Adjusters

Motoeized
Stages

Light Sources &
Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others



Application Systems

Optics & Optical Coatings ■

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources & Laser Safety

Index

Guide

Mirrors

Beamsplitters

Polarizers

Lenses

Multi-Element Optics

Filters

Prisms

Substrates/Windows

Optical Data

Maintenance

Selection Guide

Half Mirror Cube

Half Mirror Plate

Application Note

Beamsplitters

Harmonic Separator

Beam Samplers

Others