

## 1229C LaserSource Cable, 20 A, 2 m, Pigtailed



The 1229C is a 20A-rated cable for use with the 4320 and other high current instruments. The instrument end as the standard 9W4 connector, while the mount end of the cable has unterminated tinned wires for connections to custom mounts.

A version of this cable, the [1228C LaserSource Cable](#), is available which 9W4 connectors on both ends of the cable.

Only laser anode and laser cathode connections are made with this cable.

### Specifications

Instrument connection	9W4, male
Mount connection	Bare wires, tinned
Signal connections	LD anode / cathode in twisted pair PD anode / cathode in twisted pair
Current rating (A)	20
Length (m)	2
Hood construction	Metal