

1260B TECSorce Cable, 5 A, 2 m



Our standard TECSorce interconnect cable. This cable is good for applications up to 5A and includes wires for the fan connection. The cable is terminated with DB-15 connectors on both ends for quick connection to any Arroyo Instruments mount.

A version of this cable, the [1261B TECSorce Cable](#), is available which features bare tinned wires on the fixture end.

If you are using a 5300 Series TECSorce and 4-wire RTD sense measurements, or using TEC currents up to 10A, use the [1262B TECSorce Cable](#), which supports up to 10A and includes the wiring the remote sensor signals.

Specifications

| | |
|-----------------------|---------------------|
| Instrument connection | DB-15, male |
| Mount connection | DB-15, female |
| Signal connections | TE Sensor Fan |
| Current rating (A) | 5 |
| Length (m) | 2 |
| Hood construction | Metal |