

1263B TECSource Cable, 10 A, Pigtailed, 2 m



Similar to the [1261B TECSource Cable](#), but adds wiring for both the 4-wire RTD connection as well as operation to 10A. The cable is terminated with a DB-15 connector on the instrument end and bare tinned wires on the mount end for connections to non-Arroyo mounts.

A version of this cable, the [1262B TECSource Cable](#), is available which features DB-15 connectors on both ends of the cable.

DB-15 Male Connector

Pin	Color	Function / Description
1	Red	TE+
2	Red/White	TE+
3	Black	TE-
4	Black/White	TE-
5	Shield	Earth Ground
6	Shield	Earth Ground
7	Green	Sensor+
8	White	Sensor-
9	Red/Yellow	TE+
10	Black/Yellow	TE-
11	Yellow	Fan+
12	Blue	Fan-
13	n/c	No connection
14	Orange	Rem Sensor+
15	Brown	Rem Sensor-
Shell	Shield	Earth Ground

Specifications

Instrument connection	DB-15, male
Mount connection	Bare wires, tinned
Signal connections	TE Sensor Remote Sensor Fan
Current rating (A)	10
Length (m)	2
Hood construction	Metal