

284-03-150 TECMount, M2.5



The **284 TECMount** provides a flexible heating and cooling platform designed to meet demanding temperature control requirements. With a large 3" diameter cold plates with a breadboard-style mounting system, the 284 is easy to integrate into a broad range of applications.

High Temperature Option



Caution
Hot Surface!

The 284-03-150 high temperature configuration, allows for operation up to 150°C, while retaining the temperature range and thermal capacity of the standard mount. Because thermistors do not perform well at high temperatures, the sensor is replaced with a high accuracy Pt 100 RTD sensor.

- **30W thermal capacity**
 - Up to +150°C Operation
 - Active TEC cooling
 - Large cold plate w/ M2.5 holes
- Customized mounting plates available*

Flexible Cold Plate

The 284-03 cold plate features a grid of M2.5 holes plus a mounting pattern for STAR LEDs.

Customizing the Cold Plate

When our standard breadboard plate just doesn't fit, the cold plate on the 284 can be custom machined to fit the exact hole pattern of your application. We can put mounting holes just where you need them so your device mounts directly to the plate, without the need for adapters or modification.

M-Series Mounting System

An optional accessory, the **284-MKIT** is a series of mounting accessories to

accommodate the integration of the 284 mounts onto your optical bread board or other mechanical system.



Flexible Temperature Feedback

Many laser and LED devices feature an internal temperature sensor. The 284 has an external sensor input, allowing for easy switching between the plate sensor and the sensor located in the laser or LED.

Specifications

Key Specs

Thermal Capacity	30W at 25°C
Cold plate	3" diameter solid copper, hard nickel plated Breadboard, M2.5 holes

Temperature Control

Temperature Range (°C)	+15 to +150
Sensor Type	Pt 100Ω RTD
Recommended Controller	5300-08-24
TE Current, Max (A)	7.4
TE Voltage, Max (V)	14.7 (at 25°C)

General

Size (H x W x D) [in(mm)]	4.75 (120.7) x 4.4 (111.8) x 4.4 (111.8)
TEC Connector	DB15, male Phoenix 6-pin