

## Micro 9.0 9mm TO-Can LaserMount



- Inexpensive
- 9mm TO-can lasers
- Post-mountable
- Pairs with 1224 cables

### Inexpensive Mount

The Micro 9.0 was design to provide an inexpensive solution for holding and mounting 9.0mm TO-can lasers.

The Micro mounts are not intended to provide heat sinking (although they will draw a small amount of heat away), but rather a solution for post mounting a TO-can on an optical table or other test setup.

### Cables

The 1224 series of cables work with the Micro mounts to provide the electrical connections to the laser. The cables feature the standard DB-9 connection to the laser driver for easy connection to any LaserSource laser driver. The Micro mount, coupled with the appropriate 1224 cable provides a quick and inexpensive solution to your TO-can control. The third thumbnail to the left shows the mount with post (not provided), and cable (optional).

### Looking for a 5.6mm holder?

If you are looking for the same feature as the Micro 9.0, but for 5.6mm TO-can lasers, see the [Micro 5.6 LaserMount](#).

## Specifications

Laser Package Supported	9mm TO-Can
Size (H x W x D) [in(mm)]	1.0 (25.4) x 1.0 (25.4) x 0.4 (10.2)
Mounting hole	8-32 threaded hole