


## E25.M7P06.1002 (E25 connecting module)



Article number: E25M7P061002  
GTIN: 4050369052178  
Unit: piece  
Weight incl. packaging: 14,60 kg  
Weight excl. packaging: 14,60 kg

### Highlights

The E25 modular folding wall connecting module E25.M7P06.1002 from laservision is used for the individual expansion of the basic module E25.M7P06.1001.

Each individual module has the dimensions: 2100x1257x308mm (HBT) The intermediate module is filled with the proven plate material M7P06 as standard. This is based on a special aluminum sandwich structure with coated aluminum cover plates (black).

With the basic module it can be set up together with the respective number of intermediate modules, for example as a room divider in laboratories or laser-safe partitions around optical tables.

<b>COLOUR:</b>	black silver
<b>PROTECTION RANGE:</b>	Infrared near infrared
<b>PROTECTION CLASS / NORM:</b>	EN 12254 EN 60825-4
<b>PROPERTIES:</b>	both-sides usable freely selectable number of modules proven hinge system standard plate material M7P06