


# uvex

Visors & Face Shields

02.04.2020

## Product data sheet model 9710 514



**UVEX ARBEITSSCHUTZ GMBH, Würzburger Str. 181-189,  
D-90766 Fürth**

### Product features



- Protective face shield for single use
- one size fits all head sizes in the range 54-64 cm
- Fixation on the head with elastic band
- PET lens, 0.3 mm, inside fog-free
- Soft foam inner padding ensures high wearing comfort even when worn for longer periods
- Helps to keep the face free from dirt and liquids
- Suitable to keep the face free from contamination of body fluids such as blood
- Not suitable for protection against mechanical hazards, such as flying particles with higher energy
- Space-saving flat-fold packaging

<b>Technical details</b>	
<b>Headband</b>	
Material	Textile elastic band
Colour	white
<b>Lens</b>	
Material	PET
Thickness	0,3 mm
Coating	Inside anti-fog
Colour	clear
<b>Marking lens</b>	
Certification mark	CE
<b>Miscellaneous</b>	
Dimensions	330 x 230 mm (w x l)
Weight	50 g
Packaging	Individually packed in polybags 220 face shields per box Shipping carton: 69 x 38 x 35 cm
Standard	EN 166
PPE category	category II

**Warning:**

- This face shield protects against drops and splashes of liquids. Thus, it reduces the risk of a droplet infection. Due to its design, complete protection against viruses, bacteria or other biological agents is only given to a limited extent.
- This face shield is only to be used by medical personnel or rescue workers and only for the duration of the current health threat.
- Can be combined with face protection masks.
- The product must be disposed of after a single use!