

XYZ Axis Stainless Steel Extended Contact Translation Stages (Vertical)

Stage Size 25 × 25 mm / 40 × 40 mm / 60 × 60 mm / 65 × 65 mm


TSDS-5L

RoHS

XYZ axis Stainless Steel Linear Stages suitable for hard vacuum and clean room applications. Same design and performance as TSD stages.

Completely made of stainless steel for excellent corrosion resistance.

Both the stage and the stainless steel micrometer head are lubricated with a high grade vacuum grease (YVAC2).



Guide

- The TSDS series stages are suitable for use in clean rooms or vacuum applications.
- [Reference](#) Vacuum measurement data E205
- Also ideal for use in airtight spaces that have to avoid outgassing.

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motorized Stages

Light Sources & Laser Safety

Index

Guide

X Axis Stages

XY Axis Stages

Z Axis Stages

XZ Axis Stages

XYZ Axis Stages

Rotation Stages

Goniometer

Tilt Stages

Vacuum

Ball Bearing Guide

Crossed Roller

Dovetail

Lapping

V Groove Screw

Others

15 × 15 mm

25 × 25 mm

40 × 40 mm

60 × 60 mm

65 × 65 mm

80 × 80 mm

100 × 100 mm

120 × 120 mm

Others


Specifications

Part Number	METRIC	TSDS-255CL	TSDS-255SL	TSDS-405CL	TSDS-405SL	TSDS-605CL	TSDS-605SL	TSDS-655CL-M6	TSDS-655SL-M6
	INCH	—	—	TSDS-405CLUU	TSDS-405SLUU	—	—	TSDS-655CLUU	TSDS-655SLUU
Stage Size [mm]		25×25	25×25	40×40	40×40	60×60	60×60	65×65	65×65
Axes of Travel		XYZ axis							
Micrometer Position		Center	Side	Center	Side	Center	Side	Center	Side
Travel [mm]		XY axis ±3 Z axis ±3	XY axis ±3 Z axis ±3	XY axis ±6.5 Z axis ±6.5	XY axis ±6.5 Z axis ±6.5	XY axis ±6.5 Z axis ±6.5	XY axis ±6.5 Z axis ±6.5	XY axis ±6.5 Z axis ±6.5	XY axis ±6.5 Z axis ±6.5
Lead of Actuator [mm/rotation]		0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Micrometer Readable Resolution [mm]		0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
Guide Method		Extended contact ball bearing guide							
Primary Material		Stainless steel							
Finish		None							
Load Capacity [N]		14.7 (1.5kgf)	14.7 (1.5kgf)	49 (5kgf)	49 (5kgf)	49 (5kgf)	49 (5kgf)	49 (5kgf)	49 (5kgf)
Travel Accuracy	Straightness [μm]	1	1	1	1	1	1	1	1
Max. Moment Capacity	Pitch [N·m]	2.2	2.2	3.9	3.9	8.0	8.0	8.0	8.0
	Roll [N·m]	2.1	2.1	3.9	3.9	8.0	8.0	8.0	8.0
Moment Stiffness	Pitch [°/N·cm]	6.7	6.7	1.26	1.26	0.58	0.58	0.58	0.58
Weight [kg]		0.2	0.2	0.6	0.6	1.7	1.7	2.0	2.0




Outline Drawing

TSDS-255CL Hexagon socket head cap screw (stainless steel construction) M2x4...4 screws




TSDS-255SL Hexagon socket head cap screw (stainless steel construction) M2x4...4 screws




TSDS-405CL/405CLUU

Hexagon socket head cap screw (stainless steel construction) M3x6...4 screws, Spring washers
 Hexagon socket head cap screw (stainless steel construction) 6-32UNCx3/16...4 screws (UU)



TSDS-405SL/405SLUU


Hexagon socket head cap screw (stainless steel construction) M3x6...4 screws, Spring washers
 Hexagon socket head cap screw (stainless steel construction) 6-32UNCx3/16...4 screws (UU)




Part Number	A (mm)	B (mm)	C	D	E (mm)	F (mm)
TSDS-405CL	16	32	M3	M6 P1	3.5	6
TSDS-405CLUU	15.875	31.75	6-32UNC	1/4-20UNC	4	6.5

Part Number	A (mm)	B (mm)	C	D	E (mm)	F (mm)
TSDS-405SL	16	32	M3	M6 P1	3.5	6
TSDS-405SLUU	15.875	31.75	6-32UNC	1/4-20UNC	4	6.5

TSDS-605CL Hexagon socket head cap screw (stainless steel construction) M4x6...4 screws




TSDS-605SL Hexagon socket head cap screw (stainless steel construction) M4x6...4 screws




TSDS-655CL-M6/655CLUU

Hexagon socket head cap screw (stainless steel construction) M6x6...4 screws (M6)
 Hexagon socket head cap screw (stainless steel construction) 1/4-20UNCx1/4...4 screws (UU)



TSDS-655SL-M6/655SLUU

Hexagon socket head cap screw (stainless steel construction) M6x6...4 screws (M6)
 Hexagon socket head cap screw (stainless steel construction) 1/4-20UNCx1/4...4 screws (UU)



Part Number	A (mm)	B (mm)	C	D (mm)	E (mm)
TSDS-655CL-M6	25	50	M6 P1	6.5	11
TSDS-655CLUU	25.4	50.8	1/4-20UNC	7	12

Part Number	A (mm)	B (mm)	C	D (mm)	E (mm)
TSDS-655SL-M6	25	50	M6 P1	6.5	11
TSDS-655SLUU	25.4	50.8	1/4-20UNC	7	12

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

MotORIZED Stages

Light Sources & Laser Safety

Index

Guide

X Axis Stages

XY Axis Stages

Z Axis Stages

XZ Axis Stages

XYZ Axis Stages

Rotation Stages

Goniometer

Tilt Stages

Vacuum

Ball Bearing Guide

Crossed Roller

Dovetail

Lapping

V Groove Screw

Others

15 × 15 mm

25 × 25 mm

40 × 40 mm

60 × 60 mm

65 × 65 mm

80 × 80 mm

100 × 100 mm

120 × 120 mm

Others